

SU-600MD

Stand specifically designed for use with Sony LCD flat panel medical monitors



Overview

The SU-600MD monitor stand is an optional accessory for use in healthcare environments. It's designed specifically for use with Sony LCD flat panel medical monitors.

The stylish monitor stand can be fixed to a horizontal surface or desk, or to a compatible medical device or system.

Monitor angle can be easily adjusted vertically over a -10 degree (downward) to +20 degree (upward) range by holding the top and bottom of the monitor and tilting.

The SU-600MD conforms to Video Electronics Standards Association (VESA) standards.

Technical Specifications

| General | |
|--|--|
| ● Storage/Transport Temperature | -20 °C to +60 °C |
| ● Storage/Transport Humidity | 20% to 90% |
| ● Operating/Storage/Transport Pressure | 700 hPa to 1060 hPa |
| ● Dimensions (W x H x D) *1 | 350 X 366.5 X 319mm |
| ● Description | Monitor Stand |
| ● Connector | VESA |
| ● Contents | Base unit (1) Neck unit (1) M5 x 20 Screws (4) M4 x 14 Screws (4) |
| ● Mass | Approx. 2.9kg (6 lb 6.3oz) |
| ● Supplied Accessories | Instructions for Use (1) Service Contact List (1) |
| Notes | |
| ● *1 | The values for dimensions are approximate. |

Compatible Products

Surgical Monitors



LMD-2435MD
24-inch LCD Full HD medical monitor



LMD-2735MD
27-inch LCD Full HD medical monitor



LMD-X310MD
31-inch 4K 2D LCD medical monitor